

## Background

- Flavoring chemicals are added in differing concentrations to numerous flavored e-cigarette refill solutions
- Previous qualitative studies have identified commonly used flavoring chemicals among various flavors, including those with potential inhalation concern.
- 20 of these popular flavoring chemicals have been selected for quantitative identification among e-cigarette liquids.
- This study aimed to develop a validated quantitative method to measure the concentration of these chemicals (listed in Table 1).
- This method was used to measure concentrations of these 20 chemicals among 215 commercial e-cigarette refill solutions purchased in the United States.

## Instrumentation

- Analysis was performed on an Agilent 7250 gas chromatograph quadrupole/time-of-flight mass spectrometer (GC/Q-TOF)
- Column:** DB-624 UI 30 m x 0.25 mm x 1.4  $\mu$ m
- Injection Volume:** 1  $\mu$ L
- Injection Port Liner:** Universal with glass wool, 4 mm.
- Injection:** Split 20:1
- Injection Port Temp:** 320°C
- Oven:** 60°C – hold 1 min then 30°C/min to 225°C – hold 4 min. Post run at 280°C for 1 min.
- Equilibrium time:** 1.5 min
- MS Transfer line:** 250°C
- Total Run time:** 10.5 min



Figure 1. Agilent 7250 GC/Q-TOF

## Chromatogram/Calibration

Figure 1 – Total Ion Chromatogram of Targeted Flavoring Chemicals and Internal Standards





Figure 2 – Calibration Curves for Twenty Flavoring Chemicals



Circles (black) indicate calibration standards within the instrument linear range and are included in the calculation of the curve and corresponding  $r^2$  value. Triangles (blue) indicate quality control standards within the instrument linear range.

## Method Validation

Table 1 – Method and Calibration Parameters

Compound	RT (min)	Quant Ion ( $\pm 10$ ppm)	$r^2$ Value <sup>1</sup>	RSE (%)	Calibration Range <sup>2</sup> (mg/mL)	LOQ (mg/mL)
(+)Pulegone	6.754	81.0706	0.994	6.4	0.02-10.0	0.02
2,3,5-Trimethylpyrazine	5.203	122.0834	0.993	4.9	0.02-10.0	0.02
Acetoin	3.385	88.0530	0.996	8.0	0.31-10.0	0.31
Benzaldehyde	5.144	106.0409	0.993	8.5	0.02-10.0	0.02
Benzyl Alcohol	5.656	108.0566	0.994	8.6	0.02-5.0	0.02
Butanoic Acid	3.931	60.0211	0.998	8.5	0.16-10.0	0.63
DL-Limonene	5.289	93.0705	0.994	9.5	0.08-10.0	0.08
Ethyl Maltol	6.596	140.0471	0.991	9.7	0.08-10.0	0.08
Ethyl Salicylate	6.886	120.0204	0.993	8.7	0.02-10.0	0.02
Ethyl Vanillin	8.574	137.0244	0.995	6.8	0.04-10.0	0.31
Eucalyptol	5.380	139.1129	0.993	8.1	0.04-10.0	0.04
Eugenol	7.459	164.0830	0.995	8.3	0.04-10.0	0.04
Furanol	5.737	128.0478	0.991	7.8	0.02-5.0	0.02
Isovanillin	8.749	152.0475	0.996	8.1	0.04-10.0	0.31
Maltol	6.131	126.0310	0.992	7.5	0.16-10.0	0.63
L-Menthol	6.336	123.1180	0.993	7.0	0.08-10.0	0.08
Methyl Salicylate	6.491	120.0206	0.993	5.6	0.08-10.0	0.08
trans-Cinnamaldehyde	7.071	131.0499	0.993	8.9	0.02-5.0	0.16
Triacetin	7.229	103.0398	0.998	7.2	0.04-10.0	0.08
Vanillin	8.132	152.0474	0.995	6.6	0.31-10.0	0.63

<sup>1</sup>Values quantitated using Quadratic 1/x<sup>2</sup> calibration model (Internal Standard Method)  
<sup>2</sup>Calibration levels prepared by serial dilution from 10mg/mL in 50:50 PG:VG

Table 2 – Method Validation Results

Compound	Within Batch <sup>1</sup>	Between Batch <sup>2</sup>	Average Bias (%)
	Average CV (%)	Average CV (%)	
(+)Pulegone	3.8	4.6	-1.3
2,3,5-Trimethylpyrazine	2.4	4.6	1.1
Acetoin	4.5	5.8	1.2
Benzaldehyde	3.6	5.0	-4.0
Benzyl Alcohol	4.2	8.0	1.4
Butanoic Acid	2.3	3.6	-1.4
DL-Limonene	3.6	6.6	-0.2
Ethyl Maltol	3.0	4.9	2.7
Ethyl Salicylate	2.9	3.7	0.7
Ethyl Vanillin	2.8	4.0	2.0
Eucalyptol	3.3	4.9	-1.0
Eugenol	2.6	3.4	0.7
Furanol	6.2	6.7	3.8
Isovanillin	4.2	4.7	1.3
L-Menthol	4.4	5.5	-0.6
Maltol	2.2	3.9	3.5
Methyl Salicylate	2.2	3.7	1.4
trans-Cinnamaldehyde	4.9	6.6	-4.6
Triacetin	3.4	4.2	-2.6
Vanillin	2.1	6.1	1.3

<sup>1</sup>CV reported as the average of individual within batch variation for standards prepared at concentrations of 0.04, 0.10, 0.88, 1.75, 3.50 and 7.00mg/mL (each analyzed in triplicate).  
<sup>2</sup>CV reported as the average total variation (between-batch) across 7 batches for standards prepared at concentrations of 0.04, 0.10, 0.88, 1.75, 3.50 and 7.00mg/mL (each analyzed in triplicate per batch).

## Sample Preparation

30 $\mu$ L of Calibration, QC and refill solutions


30 $\mu$ L of Internal Standard (1mg/mL)

Analyze via GC/Q-TOF

3mL of Methanol Vortex 5 min

## Flavorings in E-Cigarette Refill Solutions

Figure 3 – Flavor Category Frequency Among E-Cigarette Refill Solutions



n = 215  
Alcohol and Other Flavors account for <1% of the total e-cigarette liquids selected, which are not included here

Figure 4 – Percentage of Detections Among Flavoring Chemicals in E-Cigarette Liquids





Figure 5 – Average Concentration of Flavoring Chemicals Identified in E-Cigarette Liquids



## Conclusion

- Average coefficient of variation (CV) and bias among method validation samples were within 20%, when compared within and between batches.
- This developed method is precise and reliable when analyzed across multiple days.
- 20 common chemical additives to e-cigarette refill solutions can be accurately measured simultaneously when using this method.
- The large dynamic working range allows for a variety of concentrations to be rapidly measured in consumer products.

## Disclosures

Dr. Goniewicz reports grants from Pfizer Inc. and served as an advisory board member to Johnson & Johnson, manufacturers of smoking cessation drugs, outside the submitted work.

Research reported here was supported by the National Cancer Institute of the National Institutes of Health (NIH) and the Food and Drug Administration (FDA) Center for Tobacco Products under Award Number U54CA228110. The content is solely the responsibility of the authors and does not necessarily represent the official views of the NIH or the FDA.